

## **Oto Perfect Aqua Acida**

### **Ear spray solution, 30ml**

Oto Perfect Aqua Acida ear spray solution is a medical device which contains acetic acid, isopropyl alcohol and d-panthenol. Is it used for prevention and treatment of uncomplicated infections of the outer ear canal and during water retention in the outer ear canal.

It exerts its effect by lowering the pH value in the outer ear aisle (thus creating unfavorable conditions for the growth of microorganisms), as well as due to its bactericidal and antimycotic action.

Also, Oto Perfect Aqua Acida ear spray solution is an excellent choice in water retention in your ear because it contains isopropyl alcohol, which effectively removes excess residual water that damages the skin barrier. D-panthenol in the formulation regenerates and protects the skin of the outer ear canal.

### **The use of Oto Perfect Aqua Acida ear spray solution is recommended:**

- For the treatment of mild and moderate infections of the external ear canal of bacterial and fungal origin (otitis externa - "swimming ear", otomycosis)
- As preventive measure for swimmers, divers and people prone to external ear canal infections
- In water retention in the external ear canal that occurs after swimming, diving, showering, staying in humid and warm climates, inadequate use of water-based ear cleaners and even sweating

### **How to use Oto Perfect Aqua Acida?**

Before using for the first time, press the sprayer several times until a steady dose of spray begins to emerge from the bottle. Tilt the side head and inject a single dose of Oto Perfect Aqua Acida ear spray solution, making sure that the applicator is only slightly pressed against the earlobe. Allow to act for 5 minutes, keeping head slightly tilted. Then straighten your head, let the solution drain and gently wipe your ear. In children, gently pull the earlobe backwards and then inject the product.

For infections of the external ear canal ("swimming ear") or the presence of water in the ear, apply a single dose of spray solution (in adults as needed) 3 times a day.

To prevent infections and water retention in the external ear canal, apply a single dose of spray solution to the ear canal 2-3 times a week.

### **Advantages of applying Oto Perfect aqua acida ear spray solution:**

- ✓ Thanks to the spray form, the product can be applied in a standing or sitting position, which allows comfortable application
- ✓ The spray is evenly dosed into the auditory canal
- ✓ Each dose provides the same exact amount of preparation
- ✓ No rinsing of the ear after administration is required

### **When Oto Perfect aqua acida ear spray solution should not be used:**

- In case of perforation or lack of eardrum
- In case of ear canal injury
- In case of hypersensitivity to any of the ingredients of the product
- In children under 3 years of age

**Warning for the safe use of Oto Perfect aqua acida ear spray solution:**

Oto Perfect aqua acida ear spray solution is not a medication. In case of ear pain, seek medical advice.

Oto Perfect aqua acida ear spray solution can cause mild irritation to the outer ear canal in some people. In the event of an allergic reaction to this product, discontinue use and seek medical advice.

**Advice:**

Mechanical cleaning of the ear canal (such as cotton swabs) is not recommended as it can lead to micro-injuries that are easily infected by contact with contaminated water or dust.

**How to keep Oto Perfect aqua acida ear spray solution:**

Keep Oto Perfect aqua acida ear spray solution in original packaging, in a cool and dry place, on the temperature up to 25°C. Keep the bottle in upright position, tightly closed, out of the reach of children.

**Shelf life: 2 years.**

Shelf life after first opening: 30 days.

**Packaging:** Folded cardboard box with plastic spray bottle containing 30 ml solution (approx. 220 doses) and instruction manual for medical device.

**Medical device Oto Perfect aqua clean ear spray solution contains:**

**Active substances:** 1 ml of spray solution contains 433.97 mg of isopropyl alcohol, 22.83 mg of acetic acid, 4.6 mg of D-panthenol. Solution: purified water. It contains preservative: benzalkonium chloride 0.1 mg/ml

**Other products from Oto collection:** Oto Perfect olive oil ear spray solution, Oto Perfect aqua acida ear spray solution, Oto CeruClean ear spray gel.

**Date of issue of the instruction manual of the medical device:** 25.04.2019.

**Number of decision on the registration of a medical device in the Register of Medical Devices:** 515-02-00866-19-002 from 26.09.2019.

**Manufacturer:**

ESENSA d.o.o., Viline vode bb, Belgrade, Republic of Serbia

**Distributor:**

FARMALOGIST d.o.o. Beograd, Mirijeovski bulevar 3, phone no: +381 11 33 15 090