

QUICK KHAN Massager Tiger cream gel

200 ml

NO PARABEN

FOR A LIGHT STEP

Ingredients (INCI): Aqua, Camphor, Menthol, Isopropyl Myristate, Eucalyptus Globulus Leaf Oil, Glycerin, Propylene Glycol, Methyl Salicylate, Polysorbate 20, Triethanolamine, Carbomer, Phenoxyethanol, Eugenia Caryophyllus Leaf Oil, Ethylhexylglycerin.

Khan Massage Tiger Cream Gel is formulated on the basis of ingredients that are traditionally used as an aid in degenerative and inflammatory joint and musculoskeletal disorders. It is recommended for rheumatic disorders, lumbago, neuralgia and stiffness. Carefully formulated substrates allow the release of active principles deeply into subcutaneous tissue, which decreases the pain shortly after application. Due to good absorption through the skin and high concentration on the target place, the Khan Massage Tiger Cream Gel calms the inflammatory process, enhances blood circulation, providing a pleasant warmth and relaxation of the muscles.

Because of the presence of essential oils of eucalyptus and clove, as well as menthol and camphor, it helps in all conditions caused by bad circulation. By mild rubbing of the cream gel on the place of painful and tensioned muscles, it leads to relaxation, and can also be used in the case of twist, sprain or strain caused by sport or recreational activities. The camphor has a local, natural anesthetic effect, which contributes to the pain decreasing and the overall effectiveness of the preparation in the problems encountered with the island and bruises.

In its composition, Khan Massage Tiger Cream Gel also contains methyl salicylate, which is used to locally relieve pain in rheumatic disorders. Applied to the skin, causes vasodilation and skin hyperemia (increased blood flow), followed by a pleasant warmth. By blocking transmission of painful nerve impulses from deeper muscle-bone structures, it reduces pain and inflammatory processes in the muscles and joints.

Thanks to the cream gel substrate, the product can be applied to all parts of the body, including hairy and exposed parts, leaving no greasy traces.

It helps for:

- Pain
- rheumatic discomfort
- slowed peripheral circulation
- sports injuries

How to use:

After removing the lid, top the tube from the position OFF turn into position ON. Pour the cream gel onto painful places and by circular movements with the help of massager balls on the top of the tube, rub the product up to full absorption.

Remarks:

The preparation is intended only for external use. Avoid contact with eyes, mucous membranes, or damaged skin. In case of contact, rinse immediately with plenty of water. It is not recommended to use in people who are sensitive to any of the ingredients.

Wash hands: Store at temperatures up to 25 ° C, out of the reach of children.

Period of use: 3 years.

Packaging: 200 ml cream gel in the tub with the massager.

Manufacturer: Esensa DOO Slobodna Zona Beograd (Belgrade), Viline vode bb, Republic of Serbia

Medical Marketing: Tel: +381 11 2070 965

e-mail: marketing@esensa.co.rs