

3 tubes (5g) for adult 3 tubes (5g) for children

EN Composition:

Aqua, Sorbitol, Glycerin, Propylene Glycol, Carbomer, Triethanolamine, Propylparaben, Methylparaben

What is Mikrolax?

Mikrolax gel is a product formulated based on glycerin.

How does glycerin work?

Glycerin acts as an emollient that softens, moistens, and reduces skin itching and flaking. It provides a pleasant feeling of cooling and therefore it is used for the treatment of various skin irritations which are accompanied by unpleasant difficulties. Glycerin helps to heal wounds and other changes on the skin. It does not cause irritation and can be used even on sensitive skin. It is quickly absorbed and doesn't leave a greasy residue. Due to glycerin's ability to osmotically draw water from surrounding tissues into the lumen of the intestine, it results in faster peristalsis of the intestines and their emptying. The intestine content becomes softer, and discharging the content becomes easier.

Instructions for use:

Squeeze out the contents of the tube onto the target area. If used as an enema, this process may be , as required, repeated throughout the course of the day.

Note:

The preparation should not be used in acute hemorrhoidal infections accompanied by hemorrhage, infection, and pain.

Storage conditions:

Store at temperatures of up to 22°C, away from the reach of children.

Expiry date:

3 years from the date of manufacture.

The date of manufacture is contained in the serial number imprinted on the weld of the tube and printed on the box.

Packaging:

A box containing 3 tubes.

Manufacturer:

Esensa d.o.o., Viline vode bb, Belgrade, Republic of Serbia